PCF Elettronica S.r.l.

VOC/THC portable gas chromatographic analyser Model 9388

For both continuous source and ambient monitoring



NEW!

PCF's Model 9388 VOC/THC portable gas chromatographic analyser features two extractable memory cards, of which the first one is used to set cycle programs intended to change analytical conditions as to carry out by a suitable column specific hydrocarbons analysis, while the second card is dedicated to data acquisition that lately may be acquired and managed by a PC in laboratory. Such a versatile instrument is suggested for measuring at emissions VOC /THC as well as specific hydrocarbons by direct sampling from chimney or stack or, through a heat traced line, up to 40 m distance. The detector is a carbon atom counter. Sample is introduced into a micro flame fed by hydrogen and air, where the electrical charges generated by the oxidation of Cx to CO are proportional to the hydrocarbon content in the

sample. Actual concentration is computed out of a calibration by

reference gas mixture. The electrical charges are collected by two polarised electrodes and converted by an electrical circuit into an electronic signal.

Description

PCF's Mod. 9388 is intended for measurements of reactive hydrocarbons in ambient air and/or source gas by subtracting from the total concentration of hydrocarbons the methane fraction. As the catalytic separation is not always efficient, the separation technique is based on the chromatography. A sample pump on the back of the pneumatic circuit fills a calibrated capillary of about 0.6 ml, whose volume normalised against the atmospheric pressure in order to have high reproducibility samples is injected via a 10 port rotation valve into a chromatographic column filled with PQS or a most suitable substrate that allows separation of methane fraction from the total hydrocarbons. The separated CH4 is sent to FID detector whose response is memorised by the in built micro processor. In the repair phase of the 10 port rotation valve a second sample is subsequently introduced into the FID detector to measure the total quantity of hydrocarbons (THC); the relevant electronic signal is also memorised by electronics. By subtracting the signal attributed to methane fraction from the signal of total hydrocarbons the value of hydrocarbon less methane fraction in the sample is computed.



In built zeroing of the instrument, before each measuring cycle, guarantees an high zero drift stability.

All three electronic signals are visible on display as well as available on the analogue outputs of the instrument. An INTEL 80C195 microprocessor manages all functions relative to analytical sequences as well as data management.

The in built self diagnostic software package also controls

- The operative parameters as well as the alarms status.

TECHNICAL SPECIFICATIONS

- The instrument preconditioning phase.
- The flame switching ON automatic procedure.
- Automatic alarm in case of flame OUT or hydrogen leak conditions.

The alphanumeric display with four lines and 40 characters per line shows the three measured values as ppm and/or mg/m3 as well as any occasional system failure status. Through the firmware recorded on EEPROM located on the card inserted from the front panel, the instrument can easily be set for specific hydrocarbon, e.g. benzene, toluene, styrene etc., up to maximum of six species. The chance of replacing the program card as well as the chromatographic column makes Model 9388 the best solution for monitoring any type of polluting hydrocarbon species. Measured and status data are recorded on a second 64 K PCB card. The latter inserted in a PC allows any type of data management and memorisation

- Measuring ranges (CH4/THC/NMHC)		: 0-10 ppm : 0-100 ppm	0-100 ppm or mg/m^3	
		· 0-100 ppm	0-10,000 ppm or mg/m ³	
		(further ranges	on request)	
- Background noise		$\cdot 0.01 \text{ ppm}$	on request)	
- Lower Detectable Limit (LDL)		$\sim 1\%$ full scale		
- Zero stability (24 hours)		+ 0.5%		
- Span drift (24 hours)		1 ± 0.576 $1 \pm 10\%$		
- Measuring cycle		1.170 : from 20 up to 180 sec. max		
- Response time		: from 20 up to 180 sec. max		
- Linearity		· hetter than 1%	: hetter than 1% full scale	
Precision		$\cdot + 10/2$		
- Sample flow rate		= 600 cc/min		
Operating temperature		: 5 40 °C		
Digital display		. J = 40 C : $A \times A0$ characters I CD		
Analogue outputs for each component		$\cdot 0_{-10}$ VDC	(0.10 VDC) = (4.20 mA)	
Sorial output		: standard BS 232 9 pin connector		
Zero drift		: self compensated		
- Zero/Span control		· set from front panel and/or		
- Zero, span control		remote control		
- Services	Hydrogen	: 30 ml/min		
	Air	: 300 ml/min		
- Calibration, by traceable gas cylinder		: 3 ppm CH ₄ +	: 3 ppm CH_4 + 1 ppm C_3H_8 , air balance	
- Mounting		: transportable		
- Dimensions		: 150x216x360	: 150x216x360 mm (6"x8.5"x 14", HxWxD)	
- Weight		: 8.5 Kg		
- Standard power supply		: 220/110 Vac, 5	: 220/110 Vac, 50/60 Hz, 24 Vdc (from battery) 300 VA (to be specified in order)	
- Pneumatic connections		: 1/4" or 4/6 pla	: 1/4" or 4/6 plastic or SS tubings	

How to order

Code number	Description	
041 - 4001	Mod. 9388, portable VOC/THC.CH4 monitor, 220 Vac 50 Hz	
041 - 4002	Mod. 9388, portable VOC/THC/CH4 monitor, 110 Vac 60 Hz	
041 - 4031	Mini cylinder carrying basket (three cylinders)	
041 - 4013	Hydrogen cylinder, 10 1 with pressure reducer	World wide distributor
041 - 4012	Pure air cylinder, 10 1 with pressure reducer	Via E Borsa 10
041 - 4011	Calibration gas cylinder, 10 1 with pressure reducer	20052 Monza MI Italy
041 - 4021	Expendables kit	Ph +39-039-835457 Ex +39-039-2841066
041 - 4022	Spare parts kit	E-mail criccardi@sferainformatica.it

PCF Elettronica S.r.l. via F.lli Savio, 31 24040 LEVATE (BG) Ph I +39-035-594918 Fax +39-035-4549528 e-mail pcfelettronica@spm.it